



SURU INTERNATIONAL PVT. LTD.


CENTRAL VENOUS CATHETER KITS

Product Name: Central Venous Catheter Kits - SURU CENOCATH®




Item Code: 3061, 3062, 3063, 3065

Features :

- A central venous catheter is a long, thin, flexible tube used for intensive infusion and/or transfusion therapy, invasive central venous pressure measurement and for blood collection
- The catheter is made of specially formulated bio-compatible polyurethane material which provides strength during insertion and also softens at body temperature to conform to the body tissues and reduces the risk and vascular trauma
- Specially designed soft & beveled tip for smooth & easy insertion
- Sufficiently radio-opaque material of catheter with clear, definite marking facilitates correct placement of catheter tip
- Catheter with single or double or triple or quadruple lumens in various lengths
- Soft flexible J tip guide wire prevents the vessel perforation and also provides good torque to ensure firm insertion
- The product is individually packed in tray and sterilized by ethylene oxide



| | | |
|---|--|---|
|  |  |  |
| <p>Single Lumen CVC Item Code # 3061</p> | <p>Double Lumen CVC Item Code # 3062</p> | <p>Triple Lumen CVC Item Code # 3063</p> |
| | | |

| | | |
|--|---|--|
|  |  |  |
| Quad Lumen CVC Item Code # 3065 | Antimicrobial Coated Item Code # 3063C | Guide wire Syringe with straight Introducer Needle Item code # 3063 G |

Single Lumen

| Catheter Size | Lumen Size | Length (cm) | Guidewire (inch x cm) | Y Introducer Needle | Dilator | Catheter Fastener | Luer Lock Syringe | Scalpel Blade |
|---------------|------------|-------------|-----------------------|---------------------|-----------|-------------------|-------------------|---------------|
| 2Fr | 22Ga | 5 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | 5 cc | 11 |
| | | 8 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 10 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 13 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 15 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 16 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 20 | 0.018" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| 3Fr | 20Ga | 30 | 0.018" X 70cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 5 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 8 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 10 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 13 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 15 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 16 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| 4Fr | 18Ga | 20 | 0.021" X 50cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 30 | 0.021" X 70cm | 20Ga X 38 MM | 4F X 5cm | 20Ga | | |
| | | 5 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 8 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 10 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 13 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 15 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| 5Fr | 16Ga | 16 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 20 | 0.021" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 30 | 0.021" X 70cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 5 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 8 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 10 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 13 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| 6Fr | 14Ga | 15 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 6F X 7cm | 16Ga | | |
| | | 13 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 14Ga | | |
| | | 15 | 0.035" X 50cm | 18Ga X 70 MM | 7F X 10cm | 14Ga | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 7F X 10cm | 14Ga | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 7F X 10cm | 14Ga | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 7F X 10cm | 14Ga | | |

Double Lumen

| Catheter Size | Lumen Size | Length (cm) | Guidewire (inch x cm) | Y Introducer Needle | Dilator | Catheter Fastener | Luer Lock Syringe | Scalpel Blade |
|---------------|------------|-------------|-----------------------|---------------------|-----------|-------------------|-------------------|---------------|
| 4Fr | 22Ga/22Ga | 5 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | 5 cc | 11 |
| | | 8 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 10 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 13 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 15 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 16 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 20 | 0.018" X 50cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| | | 30 | 0.018" X 70cm | 20Ga X 38 MM | 5 F X 7cm | 18Ga | | |
| 5Fr | 18Ga/20Ga | 5 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 8 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 10 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 13 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 15 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 16 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 20 | 0.021" X 50cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| | | 30 | 0.021" X 70cm | 20Ga X 38 MM | 6F X 7cm | 16Ga | | |
| 7Fr | 14Ga/18Ga | 13 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 15 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| 8Fr | 14Ga/14Ga | 13 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 15 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |

Triple Lumen

| Catheter Size | Lumen Size | Length (cm) | Guidewire (inch x cm) | Y Introducer Needle | Dilator | Catheter Fastener | Luer Lock Syringe | Scalpel Blade |
|---------------|----------------|-------------|-----------------------|---------------------|-------------|-------------------|-------------------|---------------|
| 4.5Fr | 20Ga/23Ga/23Ga | 5 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | 5 cc | 11 |
| | | 8 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| | | 10 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| | | 13 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| | | 15 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| | | 16 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| | | 20 | 0.021" X 50cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| | | 30 | 0.021" X 70cm | 20Ga X 38 MM | 6Fx7cm | 18Ga | | |
| 5.5Fr | 20Ga/22Ga/22Ga | 5 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 8 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 10 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 13 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 15 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 16 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 20 | 0.021" X 50cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| | | 30 | 0.021" X 70cm | 20Ga X 38 MM | 6.5F X 10cm | 16Ga | | |
| 7Fr | 16Ga/18Ga/18Ga | 13 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 15 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 8F X 10cm | 7Fr | | |
| 8.5Fr | 16Ga/16Ga/14Ga | 13 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 15 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |

Quad Lumen

| Catheter Size | Lumen Size | Length (cm) | Guidewire (inch x cm) | Y Introducer Needle | Dilator | Catheter Fastener | Luer Lock Syringe | Scalpel Blade |
|---------------|---------------------------|-------------|-----------------------|---------------------|-----------|-------------------|-------------------|---------------|
| 8.5Fr | 16Ga / 18Ga / 18Ga / 14Ga | 5 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | 5 cc | 11 |
| | | 8 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 10 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 13 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 15 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 16 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 20 | 0.035" X 50cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |
| | | 30 | 0.035" X 70cm | 18Ga X 70 MM | 9F X 10cm | 8Fr | | |

