



## ABSORBABLE GELATIN SPONGE

Product Name: Absorbable Gelatin Sponge

Item Code: 6007

## Features :

- White, pliable, light & insoluble in water it is capable of absorbing approximately forty times its weight of well agitated oxalated whole blood or fifty times of water
- When applied to bleeding surface in amounts sufficient to cover the
  area it rapidly controlscapillary and venous bleeding, forming a stable
  adherent coagulum it can also be used in surgery of the larger vessels
  upon application of the proper surgical techniques
- When left is situ controls recurrent hemorrhages without causing foreign body reactions
- Favors the process of tissue and when implanted in tissue is completely
  absorbed in 3-6 weeks, depending upon the quantity used and the type
  of tissue in which it is placed. It leaves only a slight amount of scar tissue
- Is compatible with antibodies when desired the sponge may be immersed in aqueous solutions of penicillin or streptomycin
- Haemostatic action of the sponge is dependent upon its proteins and physical properties. In those cases where the amounts of prothrombinnoid are sufficient, the absorbable gelatin sponge should be moistened with sterile thrombin solution
- Useful in dental surgery, malignancy, abdominal surgery, orthopedic surgery, neuro surgery, genitourinary surgery, anorectal surgery, gynecologic surgery, otolaryngological surgery

## **Available Configuration:**

Size available: 80 X 50 X 10 mm, 70 X 50 X 10 mm, 20 X 20 X 7 mm

http://www.suru.com/Surgery-Wound-Care/Absorbable-Gelatin-Sponge