



SURU INTERNATIONAL PVT. LTD.

NON - ABSORBABLE SURGICAL SUTURE U.S.P. **S290DK18**
Black Braided Silk Suture

2-0 (3 Metric)
90 cm

1/2 Circle
 Taper Cutting
 18mm

STERILE EO MD
 MFG. LIC. #
 MFG/MD/2020/000238
 LOT : E25011
 : 2025.05
 : 2030.04

SURUSIL®

Manufactured by :
SURU INTERNATIONAL PVT. LTD.
 Plot # 5 & 6, Poonamia Comp., Saravali, Jawhar Rd.,
 Dahanu 401 602, MAH., INDIA • www.suru.com
 Customer Care No.: +91 22 2820 0707. Email : suru@suru.com

1007/01/06

BLACK BRAIDED SILK SUTURE

Product Name: Black Braided Silk Suture - SURUSIL®



Item Code: 1007

Features :












- Product : Black Braided Silk Suture SURUSIL® (Non-Absorbable Surgical Suture U.S.P)
- Material : 100 % Protein fibre spun by entwined thread from silk worm Bombyx Mori Larva. Treated with wax to give strength and smooth surface
- Dye Used : Hematine Black (CI : 75290)
- Colour : Black
- Coating : Bees wax

Available Configuration :

| | | | | | | | | |
|--------------|-----|-----|-----|-----|-----|----------|---|---|
| Size (U.S.P) | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 1-0 or 0 | 1 | 2 |
| Metric | 0.7 | 1 | 1.5 | 2 | 3 | 3.5 | 4 | 5 |

| NEEDLE LENGTH | NEEDLE DESCRIPTION | SUTURE LENGTH | SUTURE SIZE (USP) | | | | |
|---|----------------------------|---------------|-------------------|---------|-----------|---------|---------|
| | | | 6-0 (0.7) | 5-0 (1) | 4-0 (1.5) | 3-0 (2) | 2-0 (3) |
|  8 mm | 3/8 Circle Reverse Cutting | 38 cm | SU 5043 | SU 5095 | | | |
|  10 mm | 3/8 Circle Round Bodied | 76 cm | | SU 5095 | | | |
|  12 mm | 3/8 Circle Cutting | 70 cm | SU 5012 | | | | |
|  12 mm Double Armed | 3/8 Circle Cutting | 70 cm | SU 5012D | | | | |
|  12 mm | 3/8 Circle Round Bodied | 76 cm | | SU 5027 | | | |
|  12 mm | 3/8 Circle Reverse Cutting | 76 cm | | SU 5080 | | | |
|  13 mm | 3/8 Circle Round Bodied | 76 cm | SU 5029 | | | | |
|  16 mm | 3/8 Circle Round Bodied | 76 cm | | SU 5081 | SU 5082 | | |
|  16 mm | 3/8 Circle Cutting | 76 cm | | | SU 5001 | SU 5002 | |
|  16 mm | 1/2 Circle Round Bodied | 76 cm | | | SU 5049 | | |
|  17 mm | 1/2 Circle Taper Cutting | 76 cm | | SU 5098 | SU 5099 | | SU 5610 |
|  20 mm | Straight Cutting | 76 cm | | SU 5090 | | | |
|  20 mm | 3/8 Circle Round Bodied | 76 cm | | SU 5096 | | | |
|  20 mm | 1/2 Circle Round Bodied | 76 cm | | | SU 5086 | SU 5087 | |
|  25 mm | 1/2 Circle Round Bodied | 76 cm | | | | SU 5070 | |
|  25 mm | 1/2 Circle Taper Cutting | 76 cm | | | | | SU 5670 |
|  26 mm | 3/8 Circle Reverse Cutting | 76 cm | | | | SU 5028 | |

Abbreviation : MP : Micro Point

| NEEDLE LENGTH | NEEDLE DESCRIPTION | SUTURE LENGTH | SUTURE SIZE (USP) | | | | |
|---|----------------------------|---------------|-------------------|---------|---------|---------|---------|
| | | | 3-0 (2) | 2-0 (3) | 0 (3.5) | 1 (4) | 2 (5) |
|  30 mm | 3/8 Circle Round Bodied | 76 cm | | | SU 5017 | | |
|  30 mm | 1/2 Circle Round Bodied | 76 cm | | SU 5333 | SU 5334 | SU 5335 | |
|  30 mm Heavy | 1/2 Circle Cutting | 76 cm | | | | SU 5061 | |
|  35 mm | 1/2 Circle Round Bodied | 76 cm | SU 5083 | | | | |
|  38 mm | 1/2 Circle Reverse Cutting | 75 cm | | SU 5003 | SU 5004 | SU 5005 | SU 5006 |
|  40 mm | Straight Cutting | 76 cm | | SU 5256 | | | |
|  45 mm | 3/8 Circle Reverse Cutting | 76 cm | | SU 5036 | SU 5037 | | |
|  50 mm | Straight Cutting | 76 cm | | SU 5208 | | | |
|  50 mm | 1/2 Circle Reverse Cutting | 76 cm | | SU 5052 | | | |
|  60 mm | Straight Cutting | 76 cm | | SU 5290 | | | |
|  60 mm | 3/8 Circle Cutting | 76 cm | | | | SU 5062 | |

PRE-CUT SUTURES (NON - NEEDED)

| SUTURE LENGTH | PACKS | SUTURE SIZE (USP) | | | | | | | | |
|---------------|------------------------|-------------------|---------|-----------|---------|---------|---------|--------|--------|--------|
| | | 6-0 (0.7) | 5-0 (1) | 4-0 (1.5) | 3-0 (2) | 2-0 (3) | 0 (3.5) | 1 (4) | 2 (5) | 3 (6) |
| 2 x 75 cm | Sterile Overwrap Packs | SU 209 | SU 210 | SU 211 | SU 212 | SU 213 | SU 214 | SU 215 | SU 216 | SU 217 |
| 10 x 35 cm | Sterile Overwrap Packs | | | | SU 222 | SU 223 | SU 224 | SU 225 | | |
| 2 x 100 cm | Sterile Overwrap Packs | | | | SU 232 | SU 233 | SU 234 | | | |

SILK REELS

| SUTURE LENGTH 25 METERS | SUTURE SIZE (USP) | | | | | | | |
|-------------------------------|-------------------|---------|-----------|---------|---------|---------|--------|--------|
| | 6-0 (0.7) | 5-0 (1) | 4-0 (1.5) | 3-0 (2) | 2-0 (3) | 0 (3.5) | 1 (4) | 2 (5) |
| 6 Reels in a Box, Non-sterile | SU 819 | SU 820 | SU 821 | SU 822 | SU 823 | SU 824 | SU 825 | SU 826 |

